



C11orf53 Antibody, FITC conjugated

Product Code	CSB-PA809008LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXP5
Immunogen	Recombinant Human Uncharacterized protein C11orf53 protein (13-236AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C11orf53, C11orf53
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	C11orf53