



C11orf53 Antibody, FITC conjugated

Product Code	CSB-PA809008LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXP5
Immunogen	Recombinant Human Uncharacterized protein C11orf53 protein (13-236AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C11orf53, C11orf53
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	C11orf53
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.