



SLC35F2 Antibody, HRP conjugated

Product Code	CSB-PA809010LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXU6
Immunogen	Recombinant Human Solute carrier family 35 member F2 protein (335-374AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 35 member F2, SLC35F2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC35F2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.