



KRTAP11-1 Antibody, Biotin conjugated

Product Code	CSB-PA810266LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUC1
Immunogen	Recombinant Human Keratin-associated protein 11-1 protein (1-186AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin-associated protein 11-1 (High sulfur keratin-associated protein 11.1), KRTAP11-1, KAP11.1 KRTAP11.1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRTAP11-1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.