



UBE2DNL Antibody

Product Code	CSB-PA810279LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWF7
Immunogen	Recombinant Human Putative ubiquitin-conjugating enzyme E2 D2-like protein (1-75AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative ubiquitin-conjugating enzyme E2 D2-like protein (Ubiquitin carrier protein D2-like) (Ubiquitin-conjugating enzyme E2D N-terminal-like) (Ubiquitin-protein ligase D2-like), UBE2DNL, UBE2D2L
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	UBE2DNL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.