



VP40 Antibody, Biotin conjugated

Product Code	CSB-PA810348HD01RCJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8JPX9
Immunogen	Recombinant Reston ebolavirus Matrix protein VP40 protein (1-331AA)
Raised In	Rabbit
Species Reactivity	Reston ebolavirus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Matrix protein VP40 (Membrane-associated protein VP40), VP40
Immunogen Species	Reston ebolavirus (strain Reston-89) (REBOV) (Reston Ebola virus)
Research Area	Others
Target Names	VP40
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.