



SLC22A9 Antibody, Biotin conjugated

Product Code	CSB-PA811607LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVM8
Immunogen	Recombinant Human Solute carrier family 22 member 9 protein (45-143AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 22 member 9 (Organic anion transporter 7), SLC22A9, hOAT4 OAT7 UST3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC22A9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.