





GLB1L2 Antibody

Product Code	CSB-PA811610LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IW92
Immunogen	Recombinant Human Beta-galactosidase-1-like protein 2 protein (401-636AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-galactosidase-1-like protein 2 (EC 3.2.1.-), GLB1L2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GLB1L2

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA811610LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA811610LA01HU at dilution of 1:100