



GLB1L2 Antibody, Biotin conjugated

Product Code	CSB-PA811610LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IW92
Immunogen	Recombinant Human Beta-galactosidase-1-like protein 2 protein (401-636AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-galactosidase-1-like protein 2 (EC 3.2.1.-), GLB1L2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GLB1L2