



MRPL41 Antibody, FITC conjugated

Product Code	CSB-PA811620LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXM3
Immunogen	Recombinant Human 39S ribosomal protein L41, mitochondrial protein (14-137AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L41, mitochondrial (L41mt) (MRP-L41) (39S ribosomal protein L27 homolog) (Bcl-2-interacting mitochondrial ribosomal protein L41) (Cell proliferation-inducing gene 3 protein) (MRP-L27 homolog) (Mitochondrial large ribosomal subunit protein mL41), MRPL41, BMRP MRPL27 RPML27
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MRPL41
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.