





HERC6 Antibody

Product Code	CSB-PA812876LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVU3
Immunogen	Recombinant Human Probable E3 ubiquitin-protein ligase HERC6 protein (764-1022AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable E3 ubiquitin-protein ligase HERC6 (EC 2.3.2.26) (HECT domain and RCC1-like domain-containing protein 6) (HECT-type E3 ubiquitin transferase HERC6), HERC6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	HERC6

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA812876LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA812876LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.