



HERC6 Antibody, HRP conjugated

Product Code	CSB-PA812876LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVU3
Immunogen	Recombinant Human Probable E3 ubiquitin-protein ligase HERC6 protein (764-1022AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable E3 ubiquitin-protein ligase HERC6 (EC 2.3.2.26) (HECT domain and RCC1-like domain-containing protein 6) (HECT-type E3 ubiquitin transferase HERC6), HERC6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	HERC6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.