




CCDC103 Antibody

Product Code	CSB-PA812877LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IW40
Immunogen	Recombinant Human Coiled-coil domain-containing protein 103 protein (1-242AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CCDC103

Image



Western Blot

Positive WB detected in: Mouse heart whole cell lysate

All lanes:CCDC103 Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 28kDa

Observed band size: 28 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.