



IZUMO1 Antibody, Biotin conjugated

Product Code	CSB-PA812903LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYV9
Immunogen	Recombinant Human Izumo sperm-egg fusion protein 1 protein (22-292AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Izumo sperm-egg fusion protein 1 (Oocyte binding/fusion factor) (OBF) (Sperm-specific protein izumo), IZUMO1
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	IZUMO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.