



# Dlk2 Antibody

<b>Product Code</b>	CSB-PA812967LA01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8K1E3
<b>Immunogen</b>	Recombinant Mouse Protein delta homolog 2 protein (27-305AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Protein delta homolog 2 (DLK-2) (Endothelial cell-specific protein S-1) (Epidermal growth factor-like protein 9) (EGF-like protein 9), Dlk2, Egfl9
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Others
<b>Target Names</b>	Dlk2

## Image



Western blot  
All lanes: Dlk2 antibody at 5µg/ml + Mouse kidney tissue  
Secondary  
Goat polyclonal to rabbit IgG at 1/10000 dilution  
Predicted band size: 41, 46, 45 kDa  
Observed band size: 41 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.