



Dlk2 Antibody

Product Code	CSB-PA812967LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8K1E3
Immunogen	Recombinant Mouse Protein delta homolog 2 protein (27-305AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein delta homolog 2 (DLK-2) (Endothelial cell-specific protein S-1) (Epidermal growth factor-like protein 9) (EGF-like protein 9), Dlk2, Egfl9
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Dlk2

Image



Western blot
All lanes: Dlk2 antibody at 5µg/ml + Mouse kidney tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 41, 46, 45 kDa
Observed band size: 41 kDa