



KRTAP8-1 Antibody, Biotin conjugated

Product Code	CSB-PA814201LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUC2
Immunogen	Recombinant Human Keratin-associated protein 8-1 protein (1-63AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin-associated protein 8-1 (High glycine-tyrosine keratin-associated protein 8.1), KRTAP8-1, KAP8.1 KRTAP8.1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KRTAP8-1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.