



MMD2 Antibody, Biotin conjugated

Product Code	CSB-PA814229LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IY49
Immunogen	Recombinant Human Monocyte to macrophage differentiation factor 2 protein (1-38AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monocyte to macrophage differentiation factor 2 (Progestin and adipoQ receptor family member 10) (Progestin and adipoQ receptor family member X), MMD2, PAQR10
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MMD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.