



TBC1D3 Antibody, Biotin conjugated

Product Code	CSB-PA814244LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZP1
Immunogen	Recombinant Human TBC1 domain family member 3 protein (1-291AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 3 (Prostate cancer gene 17 protein) (Protein TRE17-alpha) (Rab GTPase-activating protein PRC17), TBC1D3, PRC17 TBC1D3A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TBC1D3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.