



SLITRK4 Antibody, Biotin conjugated

Product Code	CSB-PA815560LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IW52
Immunogen	Recombinant Human SLIT and NTRK-like protein 4 protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SLIT and NTRK-like protein 4, SLITRK4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLITRK4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.