



NRSN1 Antibody, FITC conjugated

Product Code	CSB-PA815587LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZ57
Immunogen	Recombinant Human Neurensin-1 protein (1-65AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurensin-1 (Neuro-p24) (Vesicular membrane protein of 24 kDa) (Vesicular membrane protein p24), NRSN1, VMP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NRSN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.