



ZBTB8OS Antibody, Biotin conjugated

Product Code	CSB-PA816894LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWT0
Immunogen	Recombinant Human Protein archease protein (2-136AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein archease (Protein ZBTB8OS) (Zinc finger and BTB domain-containing opposite strand protein 8), ZBTB8OS, ARCH
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZBTB8OS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.