



NKAIN4 Antibody, FITC conjugated

Product Code	CSB-PA818240LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVV8
Immunogen	Recombinant Human Sodium/potassium-transporting ATPase subunit beta-1-interacting protein 4 protein (83-150AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium/potassium-transporting ATPase subunit beta-1-interacting protein 4 (Na(+)/K(+)-transporting ATPase subunit beta-1-interacting protein 4) (Protein FAM77A), NKAIN4, C20orf58 FAM77A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NKAIN4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.