



TRIM48 Antibody, Biotin conjugated

Product Code	CSB-PA818248LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWZ4
Immunogen	Recombinant Human Tripartite motif-containing protein 48 protein (1-208AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 48 (RING finger protein 101), TRIM48, RNF101
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRIM48
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.