



TRIM48 Antibody, Biotin conjugated

Product Code	CSB-PA818248LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWX4
Immunogen	Recombinant Human Tripartite motif-containing protein 48 protein (1-208AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 48 (RING finger protein 101), TRIM48, RNF101
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRIM48