



CDH26 Antibody

Product Code	CSB-PA818250LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXH8
Immunogen	Recombinant Human Cadherin-like protein 26 protein (1-124AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cadherin-like protein 26 (Cadherin-like protein VR20), CDH26
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CDH26
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.