



SLC38A6 Antibody, Biotin conjugated

Product Code	CSB-PA818274LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZM9
Immunogen	Recombinant Human Probable sodium-coupled neutral amino acid transporter 6 protein (132-169AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable sodium-coupled neutral amino acid transporter 6 (N-system amino acid transporter 1) (NAT-1) (Na(+)-coupled neutral amino acid transporter 6) (Solute carrier family 38 member 6), SLC38A6, NAT1 SNAT6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC38A6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.