



SLC38A6 Antibody, Biotin conjugated

Product Code	CSB-PA818274LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZM9
Immunogen	Recombinant Human Probable sodium-coupled neutral amino acid transporter 6 protein (132-169AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable sodium-coupled neutral amino acid transporter 6 (N-system amino acid transporter 1) (NAT-1) (Na(+)-coupled neutral amino acid transporter 6) (Solute carrier family 38 member 6), SLC38A6, NAT1 SNAT6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC38A6