





VKORC1L1 Antibody

Product Code	CSB-PA818673LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N0U8
Immunogen	Recombinant Human Vitamin K epoxide reductase complex subunit 1-like protein 1 protein (38-91AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vitamin K epoxide reductase complex subunit 1-like protein 1 (VKORC1-like protein 1) (EC 1.17.4.4), VKORC1L1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VKORC1L1

Image



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA818673LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human stomach tissue using CSB-PA818673LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.