



VKORC1L1 Antibody, HRP conjugated

Product Code	CSB-PA818673LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N0U8
Immunogen	Recombinant Human Vitamin K epoxide reductase complex subunit 1-like protein 1 protein (38-91AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vitamin K epoxide reductase complex subunit 1-like protein 1 (VKORC1-like protein 1) (EC 1.17.4.4), VKORC1L1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VKORC1L1