



VKORC1L1 Antibody, HRP conjugated

Product Code	CSB-PA818673LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N0U8
Immunogen	Recombinant Human Vitamin K epoxide reductase complex subunit 1-like protein 1 protein (38-91AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vitamin K epoxide reductase complex subunit 1-like protein 1 (VKORC1-like protein 1) (EC 1.17.4.4), VKORC1L1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VKORC1L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.