



CCDC155 Antibody

Product Code	CSB-PA818716LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6L0
Immunogen	Recombinant Human Protein KASH5 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein KASH5 (Coiled-coil domain-containing protein 155) (KASH domain-containing protein 5), CCDC155, KASH5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CCDC155
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.