



TTLL6 Antibody, FITC conjugated

Product Code	CSB-PA818727LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N841
Immunogen	Recombinant Human Tubulin polyglutamylase TTLL6 protein (121-310AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin polyglutamylase TTLL6 (EC 6.-.-.) (Tubulin--tyrosine ligase-like protein 6), TTLL6, TTL.6
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TTLL6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.