




NETO2 Antibody

Product Code	CSB-PA818759LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NC67
Immunogen	Recombinant Human Neuropilin and tolloid-like protein 2 protein (378-525AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuropilin and tolloid-like protein 2 (Brain-specific transmembrane protein containing 2 CUB and 1 LDL-receptor class A domains protein 2), NETO2, BTCL2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NETO2

Image



Western blot

All lanes: NETO2 antibody at 4µg/ml

Lane 1: U251 whole cell lysate


Lane 2: MDA-MB-231 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 60, 23, 59 kDa

Observed band size: 70 kDa



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA818759LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.