



NETO2 Antibody, HRP conjugated

Product Code	CSB-PA818759LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NC67
Immunogen	Recombinant Human Neuropilin and tolloid-like protein 2 protein (378-525AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuropilin and tolloid-like protein 2 (Brain-specific transmembrane protein containing 2 CUB and 1 LDL-receptor class A domains protein 2), NETO2, BTCL2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NETO2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.