



CLRN3 Antibody, Biotin conjugated

Product Code	CSB-PA818764LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NCR9
Immunogen	Recombinant Human Clarin-3 protein (39-91AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Clarin-3 (Transmembrane protein 12) (Usher syndrome type-3A-like protein 1), CLRN3, TMEM12 USH3AL1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CLRN3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.