



CLRN3 Antibody, Biotin conjugated

Product Code	CSB-PA818764LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NCR9
Immunogen	Recombinant Human Clarin-3 protein (39-91AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Clarin-3 (Transmembrane protein 12) (Usher syndrome type-3A-like protein 1), CLRN3, TMEM12 USH3AL1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CLRN3