



# MRPL30 Antibody

|                     |                                                                                                                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code        | CSB-PA819459LA01HU                                                                                                                                                                                   |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                        |
| Uniprot No.         | Q8TCC3                                                                                                                                                                                               |
| Immunogen           | Recombinant Human 39S ribosomal protein L30, mitochondrial protein (35-161AA)                                                                                                                        |
| Raised In           | Rabbit                                                                                                                                                                                               |
| Species Reactivity  | Human                                                                                                                                                                                                |
| Tested Applications | ELISA                                                                                                                                                                                                |
| Form                | Liquid                                                                                                                                                                                               |
| Conjugate           | Non-conjugated                                                                                                                                                                                       |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                                                                                                                     |
| Purification Method | >95%, Protein G purified                                                                                                                                                                             |
| Isotype             | IgG                                                                                                                                                                                                  |
| Clonality           | Polyclonal                                                                                                                                                                                           |
| Alias               | 39S ribosomal protein L30, mitochondrial (L30mt) (MRP-L30) (39S ribosomal protein L28, mitochondrial) (L28mt) (MRP-L28) (Mitochondrial large ribosomal subunit protein uL30m), MRPL30, MRPL28 RPML28 |
| Immunogen Species   | Homo sapiens (Human)                                                                                                                                                                                 |
| Research Area       | Tags & Cell Markers                                                                                                                                                                                  |
| Target Names        | MRPL30                                                                                                                                                                                               |
| Usage               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                                                                                                                          |