



LRRC15 Antibody, Biotin conjugated

Product Code	CSB-PA819484LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TF66
Immunogen	Recombinant Human Leucine-rich repeat-containing protein 15 protein (301-538AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat-containing protein 15 (Leucine-rich repeat protein induced by beta-amyloid homolog) (hLib), LRRC15, LIB
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LRRC15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.