




TGFBRAP1 Antibody

Product Code	CSB-PA819888LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUH2
Immunogen	Recombinant Human Transforming growth factor-beta receptor-associated protein 1 protein (401-600AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor-beta receptor-associated protein 1 (TGF-beta receptor-associated protein 1) (TRAP-1) (TRAP1), TGFBRAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TGFBRAP1

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA819888LA01HU at dilution of 1:100