



TGFBRAP1 Antibody, FITC conjugated

Product Code	CSB-PA819888LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUH2
Immunogen	Recombinant Human Transforming growth factor-beta receptor-associated protein 1 protein (401-600AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor-beta receptor-associated protein 1 (TGF-beta receptor-associated protein 1) (TRAP-1) (TRAP1), TGFBRAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TGFBRAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.