



LRRC29 Antibody, HRP conjugated

Product Code	CSB-PA819894LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WV35
Immunogen	Recombinant Human Leucine-rich repeat-containing protein 29 protein (1-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat-containing protein 29 (F-box and leucine-rich repeat protein 9) (F-box protein FBL9) (F-box/LRR-repeat protein 9), LRRC29, FBL9 FBXL9
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LRRC29