



LRRC29 Antibody, HRP conjugated

Product Code	CSB-PA819894LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WV35
Immunogen	Recombinant Human Leucine-rich repeat-containing protein 29 protein (1-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat-containing protein 29 (F-box and leucine-rich repeat protein 9) (F-box protein FBL9) (F-box/LRR-repeat protein 9), LRRC29, FBL9 FBXL9
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LRRC29
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.