



AFG1L Antibody, FITC conjugated

Product Code	CSB-PA819895LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WV93
Immunogen	Recombinant Human AFG1-like ATPase protein (14-313AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	AFG1-like ATPase (Lactation elevated protein 1) (EC 3.6.-.-) (Protein AFG1 homolog), AFG1L, AFG1 LACE1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.