




DYDC1 Antibody

Product Code	CSB-PA819901LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWB3
Immunogen	Recombinant Human DPY30 domain-containing protein 1 protein (1-177AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DPY30 domain-containing protein 1, DYDC1, DPY30D1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology

Image



Western blot
 All lanes: DYDC1 antibody at 4µg/ml + Hela whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 21, 22 kDa
 Observed band size: 21 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.