





# LPGAT1 Antibody

<b>Product Code</b>	CSB-PA821685LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q92604
<b>Immunogen</b>	Recombinant Human Acyl-CoA:lysophosphatidylglycerol acyltransferase 1 protein (151-341AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:3000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Acyl-CoA:lysophosphatidylglycerol acyltransferase 1 (EC 2.3.1.-), LPGAT1, FAM34A KIAA0205
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	LPGAT1

## Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA821685LA01HU at dilution of 1:100



#### Western Blot

Positive WB detected in: HEK293 whole cell lysate, HeLa whole cell lysate

All lanes: LPGAT1 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 44 kDa

Observed band size: 44 kDa