





LPGAT1 Antibody

Product Code	CSB-PA821685LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92604
Immunogen	Recombinant Human Acyl-CoA:lysophosphatidylglycerol acyltransferase 1 protein (151-341AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:3000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acyl-CoA:lysophosphatidylglycerol acyltransferase 1 (EC 2.3.1.-), LPGAT1, FAM34A KIAA0205
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LPGAT1

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA821685LA01HU at dilution of 1:100

**Western Blot**

Positive WB detected in: HEK293 whole cell lysate, HeLa whole cell lysate

All lanes: LPGAT1 antibody at 3 μ g/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 44 kDa

Observed band size: 44 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.