



RBPMs Antibody, HRP conjugated

Product Code	CSB-PA821892LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q93062
Immunogen	Recombinant Human RNA-binding protein with multiple splicing protein (1-219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein with multiple splicing (RBP-MS) (Heart and RRM expressed sequence) (Hermes), RBPMs, HERMES
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RBPMs
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.