



RBPMs Antibody, Biotin conjugated

Product Code	CSB-PA821892LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q93062
Immunogen	Recombinant Human RNA-binding protein with multiple splicing protein (1-219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein with multiple splicing (RBP-MS) (Heart and RRM expressed sequence) (Hermes), RBPMs, HERMES
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RBPMs