



ND-75 Antibody, Biotin conjugated

Product Code	CSB-PA821986LD01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q94511
Immunogen	Recombinant Drosophila melanogaster NADH-ubiquinone oxidoreductase 75 kDa subunit, mitochondrial protein (28-315AA)
Raised In	Rabbit
Species Reactivity	Drosophila melanogaster
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADH-ubiquinone oxidoreductase 75 kDa subunit, mitochondrial (EC 1.6.5.3) (EC 1.6.99.3) (Complex I-75kD) (CI-75kD), ND-75, ND75
Immunogen Species	Drosophila melanogaster (Fruit fly)
Research Area	Others
Target Names	ND-75
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.